

eIF4A1 Rabbit mAb


ABZ2351

Product Overview

Name	eIF4A1 Rabbit mAb
Accession(Primary)	P60842
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 46 kDa; Observed MW: 46 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:20-1:50

Target information (P60842)

Gene Name	EIF4A1
About Gene Symbol	EIF4A1 (Eukaryotic Translation Initiation Factor 4A1) is a metabolic enzyme involved in translational regulation. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of eIF4A1 in K562, C6, 3T3, HeLa lysates using eIF4A1 antibody.