

Rac3 Rabbit mAb


ABZ2352

Product Overview


Name	Rac3 Rabbit mAb
Accession(Primary)	P60763
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 21 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (P60763)

Gene Name	RAC3
About Gene Symbol	RAC3 (Rac Family Small GTPase 3) is a small GTPase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of RAC3 in K562, rat Brain lysates using Rac3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using RAC3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.