

PPP4C Rabbit mAb


ABZ2354

Product Overview

Name	PPP4C Rabbit mAb
Accession(Primary)	P60510
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:100; FC-1:200-1:500

Target information (P60510)

Gene Name	PPP4C
About Gene Symbol	PPP4C (Protein Phosphatase 4 Catalytic Subunit) is a protein phosphatase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PPX in K562, C6, 3T3, HeLa lysates using PPP4C antibody.