

PTEN Rabbit mAb

ABZ2355


Product Overview

Name	PTEN Rabbit mAb
Accession(Primary)	P60484
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P60484)

Gene Name PTEN

About Gene Symbol



Western blot analysis of PTEN in K562, C6, 3T3, HeLa lysates using PTEN antibody.

