

CD81 Rabbit mAb


ABZ2356

Product Overview

Name	CD81 Rabbit mAb
Accession(Primary)	P60033
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; ICC/IF-1:100; FC-1:200-1:500

Target information (P60033)

Gene Name	CD81
About Gene Symbol	



Western blot analysis of CD81 in C2C12, 3T3, MCF-7 lysates using CD81 antibody.