

YTHDF1 Rabbit mAb

ABZ2357


Product Overview

Name	YTHDF1 Rabbit mAb
Accession(Primary)	P59326
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 61 kDa; Observed MW: 61 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (P59326)

Gene Name Ythdf1

About Gene Symbol



Western blot analysis of YTHDF1 in HeLa, CHO-K1, Ramos lysates using YTHDF1 antibody.

