

SESN2/Sestrin-2 Rabbit mAb


ABZ2358

Product Overview


Name	SESN2/Sestrin-2 Rabbit mAb
Accession(Primary)	P58004
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:50-1:100; FC-1:50-1:100; IP-1:20-1:50

Target information (P58004)


Gene Name	SESN2
About Gene Symbol	



Immunocytochemistry analysis of SESN2 (green) in HeLa using SESN2 antibody, and DAPI (blue).



Western blot analysis of SESN2 in mouse spleen, C6 lysates using SESN2 antibody.



Western blot analysis of SESN2 in K562 lysates
using SESN2 antibody