


GSDMD Rabbit mAb

ABZ2359

Product Overview

Name	GSDMD Rabbit mAb
Accession(Primary)	P57764
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Recommended Dilution	WB-1:2000-1:8000

Target information (P57764)

Gene Name GSDMD**About Gene Symbol**

Western blot analysis of GSDMD in Jurkat lysates using GSDMD antibody.

