

BACE1 Rabbit mAb

ABZ2360


Product Overview

Name	BACE1 Rabbit mAb
Accession(Primary)	P56818
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 56 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:20-1:50; IP-1:20-1:50

Target information (P56818)

Gene Name Bace1

About Gene Symbol



Western blot analysis of BACE1 in rat Brain
lysates using BACE1 antibody.

