

PER3 Rabbit mAb

ABZ2362


Product Overview

Name	PER3 Rabbit mAb
Accession(Primary)	P56645
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:132kDa;Observed MW:132kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P56645)


Gene Name PER3

About Gene Symbol



Western blot analysis of PER3 in HeLa, K562 lysates using PER3 antibody.

Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of PER3 in K562 lysates
using PER3 antibody.