

## CTBP2 Rabbit mAb


ABZ2363

### Product Overview


|                             |                                                                          |
|-----------------------------|--------------------------------------------------------------------------|
| <b>Name</b>                 | CTBP2 Rabbit mAb                                                         |
| <b>Accession(Primary)</b>   | P56545                                                                   |
| <b>Host Species</b>         | Rabbit                                                                   |
| <b>Clonality</b>            | Monoclonal                                                               |
| <b>Isotype</b>              | IgG                                                                      |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                        |
| <b>Calculated MW</b>        | Calculated MW: 49 kDa; Observed MW: 49 kDa                               |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:1000-1:5000;<br>ICC/IF-1:20-1:50; FC-1:20-1:50 |

### Target information (P56545)


|                          |                                                                                                                                                                                                                                                                                                                       |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | CTBP2                                                                                                                                                                                                                                                                                                                 |
| <b>About Gene Symbol</b> | CTBP2 (C-Terminal Binding Protein 2) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |




Western blot analysis of CTBP2 in mouse lung, rat testis lysates using CTBP2 antibody.



Western blot analysis of CTBP2 in HeLa, A549, U2OS, C6 lysates using CTBP2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using CTBP2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of CTBP2 (green) in HeLa using CTBP2 antibody, and DAPI (blue).