

Bid Rabbit mAb


ABZ2364

Product Overview


Name	Bid Rabbit mAb
Accession(Primary)	P55957
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 22 kDa; Observed MW: 22 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20

Target information (P55957)


Gene Name	BID
About Gene Symbol	BID (BH3 Interacting Domain Death Agonist) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of Bid (green) in HeLa using Bid antibody, and DAPI (blue).



Western blot analysis of Bid in HeLa, A549, HL-60 lysates using Bid antibody.



Immunohistochemistry analysis of paraffin-embedded liver cancer tissue using Bid antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.