

## FOXA1 Rabbit mAb


ABZ2366

### Product Overview


|                             |                                                  |
|-----------------------------|--------------------------------------------------|
| <b>Name</b>                 | FOXA1 Rabbit mAb                                 |
| <b>Accession(Primary)</b>   | P55317                                           |
| <b>Host Species</b>         | Rabbit                                           |
| <b>Clonality</b>            | Monoclonal                                       |
| <b>Isotype</b>              | IgG                                              |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                |
| <b>Calculated MW</b>        | Calculated MW: 49 kDa; Observed MW: 49 kDa       |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC:1:20; ICC/IF-1:20; FC-1:20 |

### Target information (P55317)


|                          |                                                                                                                                                                                                                                                                                   |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | FOXA1                                                                                                                                                                                                                                                                             |
| <b>About Gene Symbol</b> | FOXA1 (Forkhead Box A1) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of FOXA1 in T47D lysates using FOXA1 antibody.



Western blot analysis of FOXA1 in C6 lysates using FOXA1 antibody.



Immunocytochemistry analysis of FOXA1 (green)  
in MCF-7 using FOXA1 antibody, and DAPI (blue)