

## VE Cadherin Rabbit mAb


ABZ2367

### Product Overview

|                             |                                                        |
|-----------------------------|--------------------------------------------------------|
| <b>Name</b>                 | VE Cadherin Rabbit mAb                                 |
| <b>Accession(Primary)</b>   | P55284                                                 |
| <b>Host Species</b>         | Rabbit                                                 |
| <b>Clonality</b>            | Monoclonal                                             |
| <b>Isotype</b>              | IgG                                                    |
| <b>Species Reactivity</b>   | Mouse                                                  |
| <b>Calculated MW</b>        | Calculated MW: 88 kDa; Observed MW: 90-140 kDa         |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; ICC/IF-1:500-1:1000;<br>IP-1:50-1:100 |

### Target information (P55284)

|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | Cdh5 |
| <b>About Gene Symbol</b> |      |



Western blot analysis of VE Cadherin in mouse lung lysates using Cadherin 5 antibody.