

R Cadherin Rabbit mAb

ABZ2368


Product Overview

Name	R Cadherin Rabbit mAb
Accession(Primary)	P55283
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 100 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:2000-1:5000


Target information (P55283)

Gene Name CDH4

About Gene Symbol



Western blot analysis of R-Cadherin/Cadherin 4 in mouse brain, U87-MG lysates using R Cadherin antibody.



Western blot analysis of RCadherin/Cadherin 4 in rat Brain, C6 lysates using RCadherin/Cadherin 4 antibody.