

ADK Rabbit mAb


ABZ2369

Product Overview


Name	ADK Rabbit mAb
Accession(Primary)	P55263
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; FC-1:200-1:500

Target information (P55263)


Gene Name	ADK
About Gene Symbol	ADK (Adenosine Kinase) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of ADK in PC-12, Raw264.7 lysates using ADK antibody.



Western blot analysis of ADK in HeLa, A549, HL-60, U251, U87-MG lysates using ADK antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using ADK antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.