

Caspase 7 Rabbit mAb


ABZ2371

Product Overview


Name	Caspase 7 Rabbit mAb
Accession(Primary)	P55210
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; FC-1:20-1:50

Target information (P55210)


Gene Name	CASP7
About Gene Symbol	



Western blot analysis of Caspase7 in Raw264.7, PC-12 lysates using Caspase 7 antibody.



Western blot analysis of Caspase7 in HeLa, A549, MCF-7, U251, HL-60 lysates using Caspase7 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Caspase7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.