

Caspase 7 Rabbit mAb

ABZ2372


Product Overview

Name	Caspase 7 Rabbit mAb
Accession(Primary)	P55210
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100


Target information (P55210)

Gene Name CASP7

About Gene Symbol



Western blot analysis of Caspase7 in A549, MCF-7, 293T lysates using Caspase 7 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using pro Caspase7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.