

PAR2 Rabbit mAb


ABZ2373

Product Overview

Name	PAR2 Rabbit mAb
Accession(Primary)	P55085
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:20-1:50; FC-1:200-1:500


Target information (P55085)

Gene Name	F2RL1
About Gene Symbol	



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of PAR2 in HepG2, MCF-7
lysates using PAR2 antibody.



Western blot analysis of PAR2 in K562 lysates
using PAR2 antibody.