

VCP Rabbit mAb


ABZ2374

Product Overview


Name	VCP Rabbit mAb
Accession(Primary)	P55072
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 89 kDa; Observed MW: 89 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:200-1:500; IP-1:10-1:100

Target information (P55072)


Gene Name	VCP
About Gene Symbol	



Immunocytochemistry analysis of VCP (green) in A549 using VCP antibody, and DAPI (blue).



Western blot analysis of VCP in MCF-7, C6, 3T3 lysates using VCP antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using VCP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.