


Aquaporin 5 Rabbit mAb

ABZ2375

Product Overview

Name	Aquaporin 5 Rabbit mAb
Accession(Primary)	P55064
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 28 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000

Target information (P55064)

Gene Name AQP5**About Gene Symbol**

Western blot analysis of Aquaporin 5 in MCF-7 lysates using Aquaporin 5 antibody.

