

Iba1 Rabbit mAb


ABZ2376

Product Overview

Name	Iba1 Rabbit mAb
Accession(Primary)	P55008
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:1000-1:4000; IHC-1:500-1:2000; ICC/IF-1:50-1:100; IP-1:20-1:50


Target information (P55008)

Gene Name	AIF1
About Gene Symbol	



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of Iba1 in THP-1 lysates using Iba1 antibody.



Western blot analysis of Iba1 in Mouse liver lysates using Iba1 antibody.