

AK2 Rabbit mAb


ABZ2377

Product Overview


Name	AK2 Rabbit mAb
Accession(Primary)	P54819
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; IP-1:10-1:100

Target information (P54819)


Gene Name	AK2
About Gene Symbol	



Immunocytochemistry analysis of AK2 (green) in Hela using AK2 antibody, and DAPI (blue).



Western blot analysis of AK2 in Hela, CHO-K1, C6, mouse Brain, Ramos lysates using AK2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using AK2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.