

## AMPK alpha 1/2 Rabbit mAb

ABZ2378

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | AMPK alpha 1/2 Rabbit mAb                  |
| <b>Accession(Primary)</b>   | P54646                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 62 kDa; Observed MW: 62 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IP-1:20-1:50              |

### Target information (P54646)

**Gene Name** PRKAA2**About Gene Symbol**

