


AMPK alpha 2 Rabbit mAb

ABZ2379

Product Overview

Name	AMPK alpha 2 Rabbit mAb
Accession(Primary)	P54646
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 62 kDa; Observed MW: 62 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (P54646)

Gene Name PRKAA2**About Gene Symbol**

Western blot analysis of AMPK alpha 2 in C6, HeLa lysates using AMPK alpha 2 antibody.

