

AMPK gamma 1 Rabbit mAb


ABZ2380

Product Overview

Name	AMPK gamma 1 Rabbit mAb
Accession(Primary)	P54619
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:1000; FC-1:20; IP-1:20

Target information (P54619)

Gene Name	PRKAG1
About Gene Symbol	PRKAG1 (Protein Kinase AMP-Activated Non-Catalytic Subunit Gamma 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and



Western blot analysis of AMPK gamma 1 in CHO-K1, C6, Jurkat lysates using AMPK gamma 1 antibody.