

AMPK gamma 1 Rabbit mAb


ABZ2381

Product Overview


Name	AMPK gamma 1 Rabbit mAb
Accession(Primary)	P54619
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:200; ICC/IF-1:50-1:100

Target information (P54619)

Gene Name	PRKAG1
About Gene Symbol	PRKAG1 (Protein Kinase AMP-Activated Non-Catalytic Subunit Gamma 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and



Western blot analysis of AMPK gamma 1 in HeLa, A549, HL-60, U2OS lysates using AMPK gamma 1 antibody.



Immunocytochemistry analysis of AMPK gamma 1 (green) in HeLa using AMPK gamma 1 antibody, and DAPI (blue)