


Ataxin 1 Rabbit mAb

ABZ2383

Product Overview

Name	Ataxin 1 Rabbit mAb
Accession(Primary)	P54253
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 87 kDa; Observed MW: 105 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:20-1:50

Target information (P54253)

Gene Name ATXN1**About Gene Symbol**

Western blot analysis of Ataxin 1 in 3T3 lysates using Ataxin 1 antibody.

