

POLG Rabbit mAb


ABZ2385

Product Overview

Name	POLG Rabbit mAb
Accession(Primary)	P54098
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 140 kDa; Observed MW: 140 kDa
Recommended Dilution	WB-1:1000-1:5000


Target information (P54098)

Gene Name	POLG
About Gene Symbol	POLG (DNA Polymerase Gamma, Catalytic Subunit) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of DNA Polymerase gamma in 293T, C6 lysates using DNA Polymerase gamma antibody.



Western blot analysis of DNA Polymerase gamma in (1) MCF-7 lysates; (2) RAW264.7 lysates using DNA Polymerase gamma antibody.