

Monocarboxylic Acid Transporter 1 Rabbit mAb

ABZ2386


Product Overview

Name	Monocarboxylic Acid Transporter 1 Rabbit mAb
Accession(Primary)	P53985
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P53985)

Gene Name SLC16A1

About Gene Symbol



Western blot analysis of Monocarboxylic acid transporter 1 in K562, HeLa lysates using Monocarboxylic Acid Transporter 1 antibody.

