

PLK1 Rabbit mAb


ABZ2388

Product Overview


Name	PLK1 Rabbit mAb
Accession(Primary)	P53350
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 68 kDa; Observed MW: 68 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:200-1:500; FC-1:50-1:100; IP-1:20-1:50

Target information (P53350)

Gene Name	PLK1
About Gene Symbol	



Western blot analysis of PLK1 in Jurkat, C6 lysates using PLK1 antibody.



Immunocytochemistry analysis of PLK1 (green) in Hela using PLK1 antibody, and DAPI (blue)