


NUP98 Rabbit mAb

ABZ2389

Product Overview

Name	NUP98 Rabbit mAb
Accession(Primary)	P52948
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 198 kDa; Observed MW: 98 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P52948)

Gene Name NUP98**About Gene Symbol**

Western blot analysis of NUP98 in HeLa, A549, HL-60 lysates using NUP98 antibody.

