

PDX1 Rabbit mAb


ABZ2390

Product Overview

Name	PDX1 Rabbit mAb
Accession(Primary)	P52945
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 31 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:20-1:100


Target information (P52945)

Gene Name	PDX1
About Gene Symbol	PDX1 (Pancreatic And Duodenal Homeobox 1) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of PDX1 in CACO-2 lysates using PDX1 antibody.

Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of PDX1 in (1) Caco2 lysates;(2) BxPC-3 lysates using PDX1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human pancreas using PDX1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.