

Hexokinase II Rabbit mAb

ABZ2391


Product Overview

Name	Hexokinase II Rabbit mAb
Accession(Primary)	P52789
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 102 kDa; Observed MW: 102 kDa
Recommended Dilution	WB-1:1000-1:2000; IHC-1:200-1:1000; IF: 1:50; FC-1:50-1:100; IP-1:20-1:50


Target information (P52789)

Gene Name	HK2
About Gene Symbol	


Rat skeletal muscle



Western blot analysis of Hexokinase II in rat skeletal muscle lysates using Hexokinase II antibody.



Western blot analysis of Hexokinase II in HeLa, CHO-K1 lysates using Hexokinase II antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Hexokinase II antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.