

VAV2 Rabbit mAb


ABZ2392

Product Overview


Name	VAV2 Rabbit mAb
Accession(Primary)	P52735
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 101 kDa; Observed MW: 101 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; FC-1:200-1:1000; IP-1:20; ICC/IF-1:50-1:100

Target information (P52735)

Gene Name	VAV2
About Gene Symbol	



Immunocytochemistry analysis of VAV2 (green) in HeLa using VAV2 antibody, and DAPI (blue).



Western blot analysis of VAV2 in HeLa, MCF-7, HepG2, U2OS, U251 lysates using VAV2 antibody.