

Eg5 Rabbit mAb


ABZ2394

Product Overview


Name	Eg5 Rabbit mAb
Accession(Primary)	P52732
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 119 kDa; Observed MW: 119 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:20-1:50; FC-1:200-1:500; IP-1:20-1:50

Target information (P52732)


Gene Name	KIF11
About Gene Symbol	KIF11 (Kinesin Family Member 11) is a kinesin motor protein involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Eg5 in HeLa, A549, HL-60, U2OS, U87-MG lysates using Eg5 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Eg5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of Eg5 (green) in
hela using Eg5 antibody, and DAPI(blue)