

MSH6 Rabbit mAb


ABZ2395

Product Overview


Name	MSH6 Rabbit mAb
Accession(Primary)	P52701
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 153 kDa; Observed MW: 163 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:200-1:500

Target information (P52701)


Gene Name	MSH6
About Gene Symbol	



Western blot analysis of MSH6 in K562, C6, 3T3, HeLa lysates using MSH6 antibody.



Immunocytochemistry analysis of MSH6(green) in HeLa using MSH6 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using MSH6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.