

MSH6 Rabbit mAb


ABZ2396

Product Overview


Name	MSH6 Rabbit mAb
Accession(Primary)	P52701
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 153 kDa; Observed MW: 163 kDa
Recommended Dilution	WB-1:500-1:1000; IHC:1:100-1:200; ICC/IF-1:20; FC-1:50-1:100; IP-1:20

Target information (P52701)

Gene Name	MSH6
About Gene Symbol	



Western blot analysis of MSH6 in A549, MCF-7, U251, C2C12 lysates using MSH6 antibody.



Immunocytochemistry analysis of MSH6 (green) in hela using MSH6 antibody, and DAPI (blue)