

hnRNP M Rabbit mAb


ABZ2398

Product Overview


Name	hnRNP M Rabbit mAb
Accession(Primary)	P52272
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 78 kDa; Observed MW: 78 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000; ICC/IF-1:100-1:200; FC-1:50-1:100

Target information (P52272)


Gene Name	HNRNPM
About Gene Symbol	HNRNPM (Heterogeneous Nuclear Ribonucleoprotein M) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of hnRNP M (green) in 293 using hnRNP M antibody, and DAPI (blue).



Western blot analysis of hnRNP M1M4 in K562, C6, 3T3, HeLa lysates using hnRNP M1M4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using hnRNP M1M4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.