

PGD Rabbit mAb


ABZ2399

Product Overview


Name	PGD Rabbit mAb
Accession(Primary)	P52209
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 53 kDa; Observed MW: 49 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:200-1:500; FC-1:100-1:1000

Target information (P52209)

Gene Name	PGD
About Gene Symbol	PGD (Phosphogluconate Dehydrogenase) is a metabolic dehydrogenase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PGD in C6, 3T3, HeLa lysates using PGD antibody.



Immunocytochemistry analysis of PGD(green) in HeLa using PGD antibody, and DAPI(blue)