

Cyclin H Rabbit mAb


ABZ2401

Product Overview


Name	Cyclin H Rabbit mAb
Accession(Primary)	P51946
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:50-1:100; FC-1:20

Target information (P51946)


Gene Name	CCNH
About Gene Symbol	



Immunocytochemistry analysis of Cyclin H (green) in HeLa using Cyclin H antibody, and DAPI (blue).



Western blot analysis of Cyclin H/p34 in 3T3 lysates using Cyclin H/p34 antibody.



Western blot analysis of Cyclin H/p34 in U251, HL-60, Lncap lysates using Cyclin H/p34 antibody