

Cyclin A2 Rabbit mAb


ABZ2402

Product Overview


Name	Cyclin A2 Rabbit mAb
Accession(Primary)	P51943
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:500; IF-1:200-1:500

Target information (P51943)

Gene Name	Ccna2
About Gene Symbol	



Western blot analysis of Cyclin A2 in Raw264.7, 3T3 lysates using Cyclin A2 antibody.



Western blot analysis of Cyclin A2 in Jurkat, C6, A431 lysates using Cyclin A2 antibody.