


RSK2 Rabbit mAb

ABZ2404


Product Overview

Name	RSK2 Rabbit mAb
Accession(Primary)	P51812
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 84 kDa; Observed MW: 84 kDa
Recommended Dilution	WB-1:1000-1:2000; IHC:1:20; FC-1:20; IP-1:20

Target information (P51812)

Gene Name RPS6KA3**About Gene Symbol**

Western blot analysis of Rsk 2/MAPKAP Kinase 1b in CHO-K1, Ramos lysates using RSK2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Rsk 2/MAPKAP Kinase 1b antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.