

STAT5b Rabbit mAb


ABZ2406

Product Overview


Name	STAT5b Rabbit mAb
Accession(Primary)	P51692
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 90 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:50-1:100; FC-1:20-1:50; IP-1:20-1:50

Target information (P51692)


Gene Name	STAT5B
About Gene Symbol	STAT5B (Signal Transducer And Activator Of Transcription 5B) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Immunocytochemistry analysis of STAT5b (green) in Hela using STAT5b antibody, and DAPI (blue).



Western blot analysis of STAT5b in Raw264.7, PC-12 lysates using STAT5b antibody.



Western blot analysis of STAT5b in HeLa, A549, HL-60 lysates using STAT5b antibody



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using STAT5b antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.