

BAP31 Rabbit mAb


ABZ2409

Product Overview


Name	BAP31 Rabbit mAb
Accession(Primary)	P51572
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 28 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:100-1:200

Target information (P51572)


Gene Name	BCAP31
About Gene Symbol	BCAP31 (B Cell Receptor Associated Protein 31) is a membrane receptor involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of BAP31 in K562, HeLa lysates using BAP31 antibody.



Immunocytochemistry analysis of BAP31(green) in HeLa using BAP31 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using BAP31 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.