

SMARCA2 / BRM Rabbit mAb


ABZ2410

Product Overview


Name	SMARCA2 / BRM Rabbit mAb
Accession(Primary)	P51531
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 181 kDa; Observed MW: 190 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:20-1:50; FC-1:200-1:500

Target information (P51531)


Gene Name	SMARCA2
About Gene Symbol	



Immunocytochemistry analysis of SMARCA2 (green) in 293 using SMARCA2 antibody, and DAPI (blue).



Western blot analysis of SMARCA2/BRM in A549, HL-60, U2OS lysates using SMARCA2/BRM antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using SMARCA2/BRM antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.