

CDKN2A/p16INK4a Rabbit mAb


ABZ2412

Product Overview


Name	CDKN2A/p16INK4a Rabbit mAb
Accession(Primary)	P51480
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 18 kDa; Observed MW: 18 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:100; FC-1:200-1:500; IP-1:20-1:50

Target information (P51480)

Gene Name	Cdkn2A
About Gene Symbol	



Western blot analysis of CDKN2A/p16INK4a in 3T3 lysates using CDKN2A antibody.



Western blot analysis of CDKN2A/p16INK4a in C6 lysates using CDKN2A/p16INK4a antibody.