

DAP1 Rabbit mAb


ABZ2413

Product Overview

Name	DAP1 Rabbit mAb
Accession(Primary)	P51397
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 11 kDa; Observed MW: 11 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200

Target information (P51397)

Gene Name	DAP
About Gene Symbol	



Western blot analysis of DAP1 in MCF-7, C6, 3T3
lysates using DAP1 antibody.