

Rab9 Rabbit mAb


ABZ2414

Product Overview


Name	Rab9 Rabbit mAb
Accession(Primary)	P51151
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; IHC-1:50-1:100

Target information (P51151)

Gene Name	RAB9A
About Gene Symbol	RAB9A (RAB9A, Member RAS Oncogene Family) is a small GTPase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Rab9 in 3T3, HeLa lysates using Rab9 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Rab9 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.