

FXR1 Rabbit mAb


ABZ2417

Product Overview


Name	FXR1 Rabbit mAb
Accession(Primary)	P51114
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70-80,82-84 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:50-1:100

Target information (P51114)


Gene Name	FXR1
About Gene Symbol	




Immunocytochemistry analysis of FXR1 (green) in Hela using FXR1 antibody, and DAPI (blue).



Western blot analysis of FXR1 in PC-12, C2C12 lysates using FXR1 antibody.



Western blot analysis of FXR1 in HeLa, MCF-7, HepG2, U2OS, U251, HL-60 lysates using FXR1 antibody



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using FXR1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.